

SHARED SPACE PACKAGE FOR AGED CARE



Soundfield Systems help achieve the best possible hearing experience for adults who are hard of hearing to;

- Improve sentence recognition ability, particularly in increased background noise
- Increase attention, interaction and participation
- Reduce Social Isolation
- Reduce activity leader vocal strain
- Can be used with a face mask
- Telecoil compatible (with Roger Neckloop)

Basic Soundfield Package for Aged Care

- 1 x Roger Touchscreen Microphone
- 1 x Roger DigiMaster 5000 V2**
- 2 x Roger Neckloop with headphones #



Full Soundfield Package for Aged Care

- 1 x Roger Touchscreen Mic
- 1 x Roger DigiMaster 5000 V2 **
- 1 x Roger Pass-Around
- 1 x Roger Media Hub
- 1 x Roger Charging Rack
- 1 x Perspex Box
- 1 x Wall Pilot
- 1 x Braille sign
- 2 x Roger Neckloop with headphones #



**** comes with tripod stand, floor stand, or wall mount—a larger sized speaker is available depending on room size**

additional Roger Neckloops can be added at any time

Recommended Applications: Dining Rooms,
Activity Spaces, Chapels, Television Rooms

DigiMaster 5000 V2 Soundfield Speaker

Delivering exceptional sound quality with the capability to produce crystal clear voices throughout a room, while reducing voice strain, distractions, background noise and feedback.

Roger Touchscreen Microphone



This easy to use microphone transmits the speaker's voice and audio feeds, whilst reducing background noise. The audio can be transmitted to the Soundfield speakers, directly to Roger receivers in hearing aids or to other portable Roger receivers such as the Roger Neckloop which links to telecoil enabled hearing aids.

- 3 intuitive directional microphones,
- 3.5mm audio input,
- Secure & private system,
- 2.4Ghz wireless transmission,

Roger Pass-Around Microphone



The Roger Pass-around Microphone is designed to enhance group discussions so speakers are heard clearly. The Roger Pass-around can be hand held or placed in the sturdy stand on a desk.

Roger Media Hub



In the modern Aged Care Facility, residents are often watching Television or listening to radio through Audiovisual sources. The Roger Multimedia Hub can be used to support both mixed (Leader and audio) and individual listening. Residents can alter volume independently of each other on the Roger Neckloop.

Perspex Box & Charging Rack



Our custom made box and lock provides secure, charging and storage for a range of devices to reduce loss and theft.

PROUDLY SUPPORTING RESIDENTS AND STAFF.

ASK US ABOUT OUR FREE TRAINING & SUPPORT SERVICES!



Hearing Loop Australia
Keeping your audience in the loop

www.hearingloop.com.au
info@hearingloop.com.au
1300 669 721