

# PORTABLE SOLUTIONS

Hearing Augmentation that provides clear, quality voice and audio feeds directly to telecoil enabled hearing aids, Roger receivers and a range of other devices.

The portable system is an expandable, comprehensive solution that is easy to use and completely wireless. Recommended for small corporate meeting and conference rooms, 1:1 meetings and small churches.

All packages are NCC D4D8 compliant, subject to signage, quantity of receivers & range requirements.



## Portable Roger Select V3 Package

The portable package includes:

- Roger Select V3,
- 2 x Roger Neckloops with headphones,
- 1 x Charging Rack with 2 neckloop adaptors;
- Portable carry case to store and protect the devices when not in use.

\* Not compatible with Soundfield Speaker & Roger Wall Pilot

## Portable Roger ON V3 Package

The portable package includes:

- Roger ON V3,
- 2 x Roger Neckloops with headphones,
- 1 x Charging Rack with 2 neckloop adaptors;
- Portable carry case to store and protect the devices when not in use.

\* Not compatible with Soundfield Speaker & Roger Wall Pilot

## Portable Roger Touchscreen V3 Package

The portable package includes:

- Roger Touchscreen Microphone V3
- 2 x Roger Neckloops with headphones,
- 1 x Charging Rack with 2 neckloop adaptors
- Portable carry case to store and protect the devices when not in use.

\* Compatible with Soundfield Speaker & Roger Wall Pilot

## Roger Select V3



This small transmitter provides a quality output to Roger receivers and hearing aids with enabled telecoils using the Roger Neckloop.

- Automatic or selectable fixed directional microphones,
- Audio Input includes Optical, 3.5mm audio input, including Bluetooth® for phone calls,
- 2.4Ghz wireless transmission,
- 10 metre range
- Unlimited Receivers for Roger Direct hearing aids.

## Roger ON V3



This small transmitter provides a quality output to Roger receivers and hearing aids with enabled telecoils using the Roger Neckloop.

- Automatic or selectable fixed directional microphones with the myRogerMic app,
- Audio Input includes Optical, 3.5mm audio input,
- Streaming audio and voice transmission during online calls,
- 2.4Ghz wireless transmission,
- 10 metre range
- Unlimited Receivers for Roger Direct hearing aids.

## Roger Touchscreen Microphone V3



This easy to use microphone transmits the speaker's voice and audio feeds, whilst reducing background noise. The audio can be transmitted to the Soundfield speakers, directly to Roger receivers in hearing aids or to the Roger Neckloop which links to telecoil enabled hearing aids.

- 3 intuitive directional microphones,
- 3.5mm audio input,
- Secure & private system,
- 2.4Ghz wireless transmission,
- 20 meter range.
- Unlimited Receivers for Roger Direct hearing aids.

## Roger Neckloop



An easy to use personal receiver, suitable for use with any hearing aid with an enabled telecoil or headphones.

- Lightweight and discrete (can be worn under clothing)
- 10 hrs battery life (comes with charger)
- Min 2 per venue (D4D8 requires 1 per 25 adults)



**Hearing Loop Australia**  
Keeping your audience in the loop

www.hearingloop.com.au  
info@hearingloop.com.au  
1300 669 721