

Auditory / Sensory Solutions

Roger Focus II



The Roger Focus II is a small, lightweight ear level receiver that sends the voice directly into the ear whilst reducing background noise.

This approach cuts out distracting noise, like nearby conversations or ambient noises allowing clearer and more comfortable understanding of speech.

Research found that children who wore Roger Focus II improved their speech understanding by an average of 53% compared to children who didn't wear the technology.

Combined with a Roger Touchscreen, ON or Select, children will have a clearer more focused opportunity to participate.

The Roger Focus II is available in 2 models—rechargeable and 312 disposable battery. The Rechargeable model is available in a choice of 10 colours, the 312 battery model is available in a choice of 2 colours.

Roger Focus II & Roger Focus II-312' Receiver Colours



Recommended Applications:

Personal communication, improved inclusion and development of social skills for Children and Teenagers with ASD, ADD, & APD

- Hear the coach clearly at footy or netball games
- Follow the instructor at dance classes
- Greater inclusion in family, community or social gatherings

Roger Touchscreen V3



Roger Touchscreen V3 is an easy to use wireless microphone. The Touchscreen has 3 intuitive directional microphones, 3.5mm audio input and an impressive 20+ metre range

Roger Select V3



Roger Select V3 is a small versatile microphone that can be worn by a presenter or placed on a table where background noise is present. When multiple conversations take place, the listener can manually select the microphone direction.

** Available in Champagne, Graphite Grey*

Roger ON V3



Roger ON V3 has the ability to automatically distinguish the direction of incoming speech. Users have the flexibility to choose the directional focus and personal settings by using the myRogerMic app. The ON V3 can be placed on a table, pointed or worn by a speaker. It has the ability to stream audio from laptops or other devices with a unique headphone feature. Great independence for teenagers. ** Available in Champagne & Graphite Grey*

Personal Auditory Package - Touchscreen

- 1 x Roger Focus II (*312 Battery or Rechargeable*)
- 1 x Roger Touchscreen Microphone V3
- * compatible with the Digimaster Soundfield Speaker*

Personal Auditory Package - Select

- 1 x Roger Focus II (*312 Battery or Rechargeable*)
- 1 x Roger Select V3

Personal Auditory Package - ON

- 1 x Roger Focus II (*312 Battery or Rechargeable*)
- 1 x Roger ON V3

For more information, scan the QR code for a link to the HLA website



Hearing Loop Australia
Keeping your audience in the loop

www.hearingloop.com.au
info@hearingloop.com.au
1300 669 721